



CD26 Polyclonal Antibody

Cat# A23518PI

-20°C/1 year

INFORMATION

Product Name	CD26 Polyclonal Antibody
Cat. No.	A23518PI
Size	100µg/50µg/20µg
Uniprot	P27487
Product type	Primary Antibody
Species Reactivity	Human
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human DPP4. AA range:350-400
Host	Rabbit
Concentration	1 mg/ml
Application	WB,IHC-p,ELISA
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year
Alias	DPP4; ADCP2; CD26; Dipeptidyl peptidase 4; ADABP; Adenosine deaminase complexing protein 2; ADCP-2; Dipeptidyl peptidase IV; DPP IV; T-cell activation antigen CD26; TP103; CD26
Image	



PRODUCT USE LIMITATION

These products are intended for research use only.