

Ovomucoid Antibody

Cat# A240001

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

INFORMATION

Product Name	Ovomucoid Antibody
Cat. No.	A240001
Size	50ug/100ug
Uniprot	P01005
Product type	Primary antibody
Species Reactivity	Chicken
Immunogen	Chicken Ovomucoid (Native Protein)
Host	Rabbit
Clonality	Polyclonal
Tested applications	ELISA
Conjugation	Non-conjugated
Purification Method	>95%, Protein G purified
Storage buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage instruction	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Alias	Ovomucoid (Allergen Gal d I) (allergen Gal d 1)

PRODUCT USE LIMITATION

These products are intended for research use only.