

SIRT3 Polyclonal Antibody

Cat# A24877PI

Storage at -20°C for 1 year

INFORMATION

Product Name	SIRT3 Polyclonal Antibody
Cat. No.	A24877PI
Size	100µg/50µg/20µg
Uniprot	Q9NTG7
Product type	Primary Antibody
Species Reactivity	Human
Immunogen	The antiserum was produced against synthesized peptide derived from human SIRT3. AA range:350-399
Host	Rb
Application	WB, ELISA
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year
Alias	SIRT3; SIR2L3; NAD-dependent protein deacetylase sirtuin-3; mitochondrial; hSIRT3; Regulatory protein SIR2 homolog 3; SIR2-like protein 3

PRODUCT USE LIMITATION

These products are intended for research use only.