

NOX-4 polyclonal antibody

Cat# A24974PI

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	NOX-4 polyclonal antibody
Cat. No.	A24974PI
Size	100µg/50µg/20µg
Uniprot	Human Q9NPH5
Product type	primary antibody
Species Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from part region of human protein
Host	Rabbit
Concentration	1 mg/ml
Clonality	polyclonal
Tested applications	WB,ELISA
Application	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	63KD
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year

PRODUCT USE LIMITATION

These products are intended for research use only.