


SARS-COV-2 Spike Monoclonal Antibody

Cat# A310002

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	SARS-COV-2 Spike Monoclonal Antibody	
Cat. No.	A310002	
Size	50ul/100ul	
Species Reactivity	SARS-COV-2	
Immunogen	Recombinant SARS-COV-2 Spike/RBD Protein (RBD) (#P0078)	
Host	Mouse / Human	
Clonity	Monoclonal	
Clone No	D003	
Tested applications	IF,ELISA	
Application	ELISA 1:5,000-10,000	
Storage buffer	0.2 µm filtered solution in PBS	
Storage instruction	Store at -20°C	
Alias	coronavirus s1,coronavirus s2,coronavirus spike,cov spike,ncov RBD,ncov s1,ncov s2,ncov spike,novel coronavirus RBD, Spike RBD	
Image		Immunofluorescence analysis of 293T cells by overexpressed SARS-CoV-2 Spike Protein (Right) or not (Left) using SARS-CoV-2 Spike Monoclonal Antibody at dilution of 1:50.

PRODUCT USE LIMITATION

These products are intended for research use only.