

SARS-COV-2 Spike Monoclonal Antibody

Cat# A310003

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	SARS-COV-2 Spike Monoclonal Antibody
Cat. No.	A310003
Size	50ul/100ul
Species Reactivity	SARS-COV-2
Immunogen	Recombinant SARS-COV-2 Spike/RBD Protein (RBD) (#P0078)
Host	Mouse / Human
Clonity	Monoclonal
Clone No	D001
Tested applications	ELISA
Application	ELISA 1:5,000-10,000
Storage buffer	0.2 µm filtered solution in PBS
Storage instruction	Store at -20°C
Alias	coronavirus s1,coronavirus s2,coronavirus spike,cov spike,ncov RBD,ncov s1,ncov s2,ncov spike, novel coronavirus RBD, novel coronavirus spike, RBD,Spike RBD

PRODUCT USE LIMITATION

These products are intended for research use only.