

SPG7 Polyclonal Antibody

Cat# A31344PI

Upon receipt, store at -20°C for 1 year

INFORMATION

Product Name	SPG7 Polyclonal Antibody
Cat. No.	A31344PI
Size	100µg/50µg/20µg
Uniprot	Q9UQ90
Product type	Primary antibody
Species Reactivity	Human
Immunogen	Synthesized peptide derived from human protein . at AA range: 71-120
Host	Rabbit
Concentration	1 mg/ml
Clonality	Polyclonal
Tested applications	WB,ELISA
Application	WB 1:500-2000 ELISA 1:5000-20000
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage instruction	-20°C/1 year

PRODUCT USE LIMITATION

These products are intended for research use only.