



Mouse Anti Human CD5 Antibody (PE-Cyanine7)

Cat# A330132

Store at 2-8°C. Avoid repeated freeze.

INFORMATION

Product Name	Mouse Anti Human CD5 Antibody (PE-Cyanine7)
Cat. No.	A330132
Size	20T/50T/100T/200T
Uniprot	P06127
Species Reactivity	Human
Host	Mouse
Clone No.	UCHT2
Tested applications	FCM
Application	The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood).
Conjugation	PE/Cyanine7
Storage buffer	PBS with 0.05% Proclin300, 1% BSA
Storage instruction	Store at 2-8°C. Avoid repeated freeze.
Alias	T-cell surface glycoprotein CD5,CD5,Lymphocyte antigen T1/Leu-1,LEU1
Image	<p style="text-align: center;">PE/Cyanine7 Excitation and Emission Spectra</p>  <p style="text-align: center;">Ex:495;565;755 nm; Em:775 nm</p>



PRODUCT USE LIMITATION

These products are intended for research use only.