

GAPDH Mouse Monoclonal Antibody

Cat# A33263PI

store at -20°C

INFORMATION

Size	100 µl
Host/Isotope	Mouse Monoclonal
Reactivity	Human,Rat,Mouse,Mk,Dg,Ch,Hamster,Rb,Pig,sheep,Insect,Yeast
Applications	WB,IF/ICC,IHC-p
Concentration	1 mg/ml
Immunogen	Synthetic Peptide of GAPDH
Specificity	The antibody detects endogenous GAPDH protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Purity	≥95% as determined by SDS-PAGE
Dilution	WB: 1:5000-20000 IHC: 1:200-300 IF 1:200
Purification	The antibody was affinity-purified from mouse ascites by affinitychromatography using specific immunogen.
Storage Stability	-20°C/1 year
Observed Band	37 KD

PRODUCT USE LIMITATION

These products are intended for research use only.