


# COX2/PTGS2 Polyclonal Antibody

Cat# A340001

Upon receipt, store at -20°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	COX2/PTGS2 Polyclonal Antibody
<b>Cat. No.</b>	A340001
<b>Size</b>	20ul/100ul
<b>Uniprot</b>	Human P35354/Mouse Q05769/Rat P35355
<b>Product type</b>	Primary antibody
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse COX2
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	0.84 mg/mL
<b>Tested applications</b>	WB
<b>Application</b>	WB 1:500-1:2000
<b>Molecular Weight</b>	74kDa
<b>Conjugation</b>	Unconjugated
<b>Purification Method</b>	Affinity purification
<b>Storage buffer</b>	PBS with 0.02% sodium azide, 1% BSA and 50% glycerol, pH7.4
<b>Storage instruction</b>	Store at -20°C . Avoid freeze / thaw cycles.
<b>Alias</b>	COX2,COX-2,COX2,Cyclooxygenase2,Cyclooxygenase 2b,Cyclooxygenase,Cyclooxygenase-2,Cyclooxygenase2,EC 1.14.99.1,fj02a10,Glucocorticoid-regulated inflammatory cyclooxygenase,Glucocorticoid-regulated inflammatory Prostaglandin G/H synthase,GRIPGHS,hCox 2,Macrophage activation-associated marker protein P71/73,OTTHUMP00000033524,PES-2,PGG/HS,PGH synthase 2,PGH2,PGHS 2,PGHS-2,PGHS2,PHS 2,PHS II,PHS2,Prostaglandin endoperoxide synthase 2 (prostaglandin G/H synthase and cyclooxygenase),Prostaglandin endoperoxide synthase 2,Prostaglandin G/H synthase 2,Prostaglandin G/H synthase 2 precursor,Prostaglandin G/H synthase and cyclooxygenase,Prostaglandin G/H synthase,Prostaglandin H2 synthase 2,prostaglandin-endoperoxide synthase 2 (prostaglandin G/H synthase and cyclooxygenase),Prostaglandin-endoperoxide synthase 2,PTGS2,ptgs2a,TIS10,TIS10 protein,unp1239,wu:fj02a10



### **PRODUCT USE LIMITATION**

These products are intended for research use only.