


Anti-SARS-CoV-2 NP Polyclonal Antibody

Cat# A340002

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-SARS-CoV-2 NP Polyclonal Antibody	
Cat. No.	A340002	
Size	50ul/100ul	
Uniprot	Human Q04206 ; Mouse Q04207	
Reactivity	SARS-COV2	
Immunogen	Recombinant SARS-COV-2 Nucleocapsid Protein (His Tag) (#P0077)	
Host	Rabbit	
Concentration	1 mg/ml	
Clonality	Polyclonal	
Tested applications	WB,ELISA	
Application	WB 1:1,000-5,000 ; ELISA 1:5,000-10,000	
Purification Method	The specific IgG was purified by Recombinant SARS-COV-2 Nucleocapsid Protein affinity purification	
Storage buffer	0.2 μm filtered solution in PBS	
Storage instruction	Store at -20°C	
Alias	coronavirus NP,coronavirus Nucleocapsid,coronavirus Nucleoprotein,cov np,ncov NP,novel coronavirus NP,novel coronavirus Nucleocapsid,novel coronavirus Nucleoprotein, NP, Nucleocapsid,Nucleoprotein Antibody	
Image		Western Blot analysis of SARS-CoV2- NP protein (#P0077) using SARS-COV-2 NP Polyclonal Antibody at dilution of 1:2000

PRODUCT USE LIMITATION

These products are intended for research use only.