

Anti-SARS-CoV-2 NP Polyclonal Antibody

Cat# A340003

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-SARS-CoV-2 NP Polyclonal Antibody
Cat. No.	A340003
Size	50ul/100ul
Reactivity	SARS-COV2
Immunogen	Recombinant SARS-COV-2 Nucleocapsid Protein (His Tag) (#P0077)
Host	Rabbit
Concentration	1 mg/ml
Clonality	Polyclonal
Tested applications	ELISA
Application	ELISA 1:5,000-10,000
Purification Method	The specific IgG was purified by Recombinant SARS-COV-2 Nucleocapsid Protein affinity purification
Storage buffer	0.2 µm filtered solution in PBS
Storage instruction	Store at -20°C
Alias	coronavirus NP, coronavirus Nucleocapsid, coronavirus Nucleoprotein, cov np, ncov NP, novel coronavirus NP, novel coronavirus Nucleocapsid, novel coronavirus Nucleoprotein, NP, Nucleocapsid, Nucleoprotein Antibody

PRODUCT USE LIMITATION

These products are intended for research use only.