



Phospho NFκB/p-NFκB p65 (Ser536) Polyclonal Antibody

Cat# A340004

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Phospho NFκB/p-NFκB p65 (Ser536) Polyclonal Antibody	
Cat. No.	A340004	
Size	20ug/50ug/100ug	
Uniprot	Q04206	
Product type	Primary antibody	
Species Reactivity	Human,Mouse,Rat	
Immunogen	A phospho specific peptide corresponding to residues surrounding S536 of human NF-κB p65	
Host	Rabbit	
Concentration	1 mg/ml	
Clonality	Polyclonal	
Tested applications	WB,IHC	
Application	WB 1:500-1:1000 IHC 1:50-1:100	
Molecular Weight	71kDa	
Purification Method	Affinity purification	
Storage buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.	
Storage instruction	-20°C/1 year	
Alias	NFKB3,p65,NF-κB p65,RELA,CMCU	
Image		<p>Western blot analysis of extracts of HeLa cells using Phospho-NF-κB p65(S536) Polyclonal Antibody at dilution of 1:1000. HeLa cells were treated by TNF-α (20ng/ml) for 30 minutes.</p>

		Immunohistochemistry of paraffin-embedded Mouse testis using Phospho-NF-kB p65(S536) Polyclonal Antibody at dilution of 1:200 (40x lens).
--	---	---

PRODUCT USE LIMITATION

These products are intended for research use only.