




NFκB p65 Polyclonal Antibody

Cat# A340005

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	NFκB p65 Polyclonal Antibody	
Cat. No.	A340005	
Size	20ug/50ug/100ug	
Uniprot	Q04206	
Product type	Primary antibody	
Species Reactivity	Human,Mouse,Rat	
Immunogen	Synthesized peptide derived from human NFκB-p65 around the non-phosphorylation site of Ser281.	
Host	Rabbit	
Concentration	1 mg/ml	
Clonality	Polyclonal	
Tested applications	WB,IHC-p,ELISA	
Application	WB 1:500-1:2000, IHC 1:100-1:300, ELISA 1:5000	
Molecular Weight	60kDa	
Purification Method	Affinity purification	
Storage buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.	
Storage instruction	-20°C/1 year	
Alias	RELA,NFKB3,Transcription factor p65,Nuclear factor NF-kappa-B p65 subunit,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3	
Image		Western Blot analysis of KB, HeLa, 293T, 3T3 cells using NFκB-p65 Polyclonal Antibody at dilution of 1:2000.

		<p>Immunohistochemistry of paraffinembedded Human stomach cancer tissue using NFκB-p65 Polyclonal Antibody at dilution of 1:200.</p>
		<p>Immunohistochemistry of paraffinembedded Mouse kidney tissue using NFκB-p65 Polyclonal Antibody at dilution of 1:200.</p>

PRODUCT USE LIMITATION

These products are intended for research use only.