


IGFBP1 Monoclonal Antibody

Cat# A34970MC

Upon receipt, store at -20°C or -80°C . Avoid repeated freeze.

INFORMATION

Product Name	IGFBP1 Monoclonal Antibody
Cat. No.	A34970MC
Size	100ug/50ug
Uniprot	P08833
Product type	Primary antibody
Species Reactivity	Human,Mouse
Immunogen	Recombinant Human IGFBP1 Protein
Host	Mouse
Tested applications	ELISA,WB
Application	WB1:500-1:5000
Clonality	monoclonal
Clone No.	1E1B2
Purification Method	>95%,Protein G purified
Storage buffer	Preservative: 15mmol/L NaN3Constituents: 0.1mol/L NaCl
Storage instruction	-20°C/1 year
Alias	AFBP antibody; Alpha pregnancy associated endometrial globulin antibody; Amniotic fluid binding protein antibody; Binding protein 25 antibody; Binding protein 26 antibody; Binding protein 28 antibody; Growth hormone independent binding protein antibody; hIGFBP 1 antibody; hIGFBP1 antibody; IBP 1 antibody; IBP-1 antibody; IBP1 antibody; IBP1_HUMAN antibody; IGF binding protein 1 antibody; IGF BP25 antibody; IGF-binding protein 1 antibody; IGFBP 1 antibody; IGFBP-1 antibody; IGFBP1 antibody; Insulin like growth factor binding protein 1 antibody; Insulin-like growth factor-binding protein 1 antibody; Placental protein 12 antibody; PP 12 antibody; PP12 antibody



PRODUCT USE LIMITATION

These products are intended for research use only.