


# CRB2 Polyclonal Antibody

Cat# A440108

Upon receipt, store at -20°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	CRB2 Polyclonal Antibody
<b>Cat. No.</b>	A440108
<b>Size</b>	100ul
<b>Uniprot</b>	Q5IJ48
<b>Product type</b>	Primary Antibody
<b>Species Reactivity</b>	Mouse, Rat ; Predicted Cross Reactive Species: Human
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human CRB2, 701-800/1285
<b>Target</b>	CRB2
<b>Host</b>	Rabbit
<b>Concentration</b>	1ug/ul
<b>Clonality</b>	Polyclonal
<b>Tested applications</b>	WB, ELISA, IHC-P, IHC-F, IF(IHC-P), IF(IHC-F), IF(ICC), ICC
<b>Application</b>	WB(1:300-1000) ELISA(1:500-1000) IHC-P(1:200-400) IHC-F(1:100-500) IF(IHC-P)(1:50-200) IF(IHC-F)(1:50-200) IF(ICC)(1:50-200) ICC(1:100-500)
<b>Conjugation</b>	Unconjugated
<b>Purification Method</b>	Purified by Protein A.
<b>Storage buffer</b>	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
<b>Storage instruction</b>	Store at -20°C for 12 months.
<b>Alias</b>	Protein crumbs homolog 2; Crumbs-like protein 2; CRB2



**PRODUCT USE LIMITATION**

These products are intended for research use only.