

Anti-Hba1C monoclonal antibody

Cat# A530110

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Anti-Hba1C monoclonal antibody
Cat. No.	A530110
Size	1 mg
Product type	Monoclonal antibody
Immunogen	Synthetic Peptide of HbA1c
Target	The antibody can detect HbA1c protein
Host	Mouse
Clonality	Monoclonal antibody
Tested applications	LETIA (Latex particle-enhanced turbidimetric immunoassay)
Conjugation	Unconjugated
Purification Method	Affinity purified
Storage instruction	Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage Buffer	10 mM PB + 150 mM NaCl + 0.1% P300, pH7.4
Concentration	Lot specific
Alias	Alpha-globin; Hemoglobin alpha chain

PRODUCT USE LIMITATION

These products are intended for research use only.