


Progesterone Receptor(PR) (6P2) mouse mAb

Cat# A630002

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Progesterone Receptor(PR) (6P2) mouse mAb	
Cat. No.	A630002	
Size	100µg/50µg	
Uniprot	P06401	
Product type	Primary antibody	
Species Reactivity	Human	
Immunogen	Synthesized peptide derived from human Progesterone Receptor(PR)	
Host	Mouse/IgG1, Kappa	
Concentration	1 mg/ml	
Clonality	Monoclonal	
Tested applications	IHC-p	
Application	IHC-p 1:100-500	
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.	
Storage instruction	-20°C/1 year	
Image		<p>Immunohistochemical analysis of paraffin-embedded Breast carcinoma-high magnification.</p> <ol style="list-style-type: none"> 1, Antibody was diluted at 1:200 (4 °C, overnight). 2, TRIS-EDTA of pH8.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 30min).

PRODUCT USE LIMITATION

These products are intended for research use only.