

Mouse Monoclonal Antibody to KI67

Cat# A630004

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Mouse Monoclonal Antibody to KI67
Cat. No.	A630004
Size	100µg/50µg
Product type	Primary antibody
Species Reactivity	Human
Immunogen	AA range : CEDLAGFKELFQTPG
Host	Mouse
Clonality	Monoclonal
Clon. No.	8D5B6
Tested applications	WB,ELISA,IHC-p
Molecular Weight	358 KD
Storage instruction	-20°C/1 year
Alias	KIA; Ki-67; MKI67

PRODUCT USE LIMITATION

These products are intended for research use only.