

ER α Polyclonal Antibody

Cat# A640001

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	ER α Polyclonal Antibody
Cat. No.	A640001
Size	100 μ g/50 μ g/20ug
Uniprot	Human P03372 / Mouse P19785 / Rat P06211
Product type	Primary antibody
Species Reactivity	Human,Mouse,Rat
Immunogen	The antiserum was produced against synthesized peptide derived from human Estrogen Receptor-alpha. AA range:71-120
Host	Rabbit
Concentration	1 mg/ml
Clonality	Polyclonal
Tested applications	WB,ELISA,IHC-p
Application	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Molecular Weight	66 KD
Storage instruction	-20°C/1 year
Alias	ESR1; ESR; NR3A1; Estrogen receptor; ER; ER-alpha; Estradiol receptor; Nuclear receptor subfamily 3 group A member 1

PRODUCT USE LIMITATION

These products are intended for research use only.