

HSP70 Polyclonal Antibody



Cat#: A20808PI

Store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	HSP70 Polyclonal Antibody
Cat. No.	A20287PI
Size	100µg/50µg/20µg
Reactivity	Human, Mouse, Rat
Gene name	HSPA1L/HSPA1A
Protein name	Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B
Human gene id	3303/3304/3005
Human swissprot	P34931/P08107
Rat gene id	24963/24472
Rat swissprot	P55063/Q07439
Immunogen	The antiserum was produced against synthesized peptide derived from human HSPA1L. AA range:221-270
Specificity	HSP70 Polyclonal Antibody detects endogenous levels of HSP70 protein..
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Rabbit
Applications	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration	1 mg/ml
Molecular Weight	70 KD

PRODUCT IMAGES



PRODUCT USE LIMITATION

These products are intended for research use only.