

Normal Rat Serum

Cat# P0068 – 10ml | P0069 – 50ml | P0070 – 100ml

Store at -20°C .

INFORMATION

Product Name	Normal Rat Serum
Cat. No.	P0068/P0069/P0070
Size	10ml/50ml/100ml
Description	The Normal Rat Serum is collected from healthy rat by aseptic collection and ensures free of bacteria, fungi, mycoplasma and virus contamination by passing three 0.1µm filters. All operations comply with the GMP (Good Manufacturing Practices for Drug) Production Standards.
Application	Cell culture, blocking step of immunofluorescence staining and immunohistochemistry et al

PRODUCT USE LIMITATION

These products are intended for research use only.