

Human IgG, Purified

Cat# P0103

Store at 2-8°C

INFORMATION

Product Name	Human IgG, Purified
Cat. No.	P0103
Size	1mg/10mg
Host	Human
Concentration	> 4.5 mg/ml (E 1% at 280 nm = 13.0)
Tested applications	Western Blot, ELISA
Conjugation	Unconjugated
Purification Method	> 95% (SDS-PAGE), Protein A
Storage buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Storage instruction	Store at 2-8°C

PRODUCT USE LIMITATION

These products are intended for research use only.