

# Human IgG, Whole Molecule

Cat# P0104

Store at 2-8°C

## INFORMATION

<b>Product Name</b>	Human IgG, Whole Molecule
<b>Cat. No.</b>	P0104
<b>Size</b>	1mg
<b>Host</b>	Human
<b>Concentration</b>	1.5 mg/ml
<b>Tested applications</b>	Western Blot, ELISA, Immunohistochemistry (IHC)
<b>Conjugation</b>	Biotin (NHS-LC)
<b>Purification Method</b>	Purified
<b>Storage buffer</b>	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free
<b>Storage instruction</b>	Store at 2-8°C

## PRODUCT USE LIMITATION

These products are intended for research use only.