

Rabbit IgG Purified, Immunoglobulin G

Cat# P0129

Upon receipt, store at 2 -8°C. Avoid repeated freeze.

INFORMATION

Product Name	Rabbit IgG Purified, Immunoglobulin G
Cat. No.	P0129
Size	10mg/50mg/100mg
Host	Rabbit
Concentration	10mg/mL
Tested applications	Western Blot, ELISA, Blocking Agent
Purification Method	Protein A
Storage buffer	10 mM Sodium Phosphate, 0.15 M Na/K Chloride, pH 7.2
Storage instruction	2-8 °C

PRODUCT USE LIMITATION

These products are intended for research use only.