

Goat Anti-Human IgG (H&L), HRP

Cat# A32206PI

Upon receipt, store at -20°C and protected from light. Avoid repeated freeze.

INFORMATION

Product Name	Goat Anti-Human IgG (H&L), HRP
Cat. No.	A32206PI
Size	100ul/1ml
Product type	Secondary Antibody
Species Reactivity	human
Host	Goat
Concentration	1-1.2mg/ml
Tested applications	ELISA, WB, IHC
Application	ELISA:1/4000-1/20000 ; WB:1/4000-1/20000 ; IHC:1/500-1/1000
Conjugation	HRP
Purification Method	Antigen affinity purification
Storage buffer	0.02M PBS, 50% glycerol, pH7.4
Storage instruction	Store at -20°C and protected from light.

PRODUCT USE LIMITATION

These products are intended for research use only.