


Donkey anti-Goat IgG (H&L), DyLight 488

Cat# A520001

Store freeze-dried powder at 2-8 °C

INFORMATION

Product Name	Donkey anti-Goat IgG (H&L) , DyLight 488
Cat. No.	A520001
Size	1 mg
Product type	Secondary antibodies
Species Reactivity	Goat
Immunogen	Purified Goat IgG, whole molecule
Host	Donkey
Concentration	1 mg/mL
Tested applications	IF, FCM, ELISA
Conjugation	DyLight 488
Purification Method	Affinity purified using solid phase Goat IgG
Storage buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free
Storage instruction	Store freeze-dried powder at 2-8 °C Product is stable for up to 4 weeks at 2-8°C after rehydration. For extended storage after rehydration, add an equal volume of glycerol and store at -20°C.
Image	

PRODUCT USE LIMITATION

These products are intended for research use only.