


# Donkey anti-Goat IgG (H&L), DyLight 594

Cat# A52002

Store freeze-dried powder at 2-8 °C

## INFORMATION

<b>Product Name</b>	Donkey anti-Goat IgG (H&L) , DyLight 594
<b>Cat. No.</b>	A520002
<b>Size</b>	1 mg
<b>Product type</b>	Secondary antibodies
<b>Species Reactivity</b>	Goat
<b>Immunogen</b>	Purified Goat IgG, whole molecule
<b>Host</b>	Donkey
<b>Concentration</b>	1 mg/mL
<b>Tested applications</b>	IF, FCM, ELISA
<b>Conjugation</b>	DyLight 594
<b>Purification Method</b>	Affinity purified using solid phase Goat IgG
<b>Storage buffer</b>	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free
<b>Storage instruction</b>	Store freeze-dried powder at 2-8 °C Product is stable for up to 4 weeks at 2-8°C after rehydration. For extended storage after rehydration, add an equal volume of glycerol and store at -20°C.
<b>Image</b>	

## PRODUCT USE LIMITATION

These products are intended for research use only.