




Myc-tag (7D9) Mouse mAb

Cat# A530104

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Myc-tag (7D9) Mouse mAb
Cat. No.	A530104
Size	50ul
Product type	Tag antibody
Immunogen	Synthetic Peptide of Myc-Tag
Target	The antibody can detect Myc fusion proteins.
Host	Mouse
Clonality	Monoclonal antibody
Clone No.	7D9-4C6-6F10
Tested applications	WB , IP , IF
Application	WB: 1/5000 IP: 1/200 IF: 1/1000
Conjugation	Unconjugated
Purification Method	Affinity purified
Storage instruction	Store at 4°C short term.Store at -20°C long term.Avoid freeze / thaw cycle.
Alias	c-Myc
Image	<p>1ug Myc fusion protein+ Primary antibody dilution at 1) 1:5000. 2) 1:10000</p>  <p>IF analysis of 293 cells transfected with a Myc-tag protein,1:2000 dilution (blue DAPI,red anti-Myc)</p> 

		<p>IP antibody use : 3ug Flag Mouse IgG1 per ml Lysate,WB 1:50001 、 untransfected 293 cell lysate2 、 transfected 293 cell lysate with C-Myc-tag fusion protein3 、 IP (transfected 293+ normal Mouse IgG+Protein G agarose)4 、 IP (transfected 293+anti- C-Myc mAb+ Protein G agarose)5 、 IP (transfected 293+Protein G)6 、 Recombinant protein (E.coli)</p>
--	---	---

PRODUCT USE LIMITATION

These products are intended for research use only.