




Flag-tag (7E1) Mouse mAb

Cat# A530105

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	Flag-tag (7E1) Mouse mAb
Cat. No.	A530105
Size	50ul
Product type	Tag antibody
Immunogen	Synthetic Peptide of Flag-Tag
Target	The antibody detects C-terminal, internal, and N-terminal Flag-tag fusion proteins.
Host	Mouse
Clonality	Monoclonal antibody
Clone No.	7E1-10D2-3H1
Tested applications	WB , IP , IF
Application	WB: 1/5000 IP: 1/200 IF: 1/2000
Conjugation	Unconjugated
Purification Method	Affinity purified
Storage instruction	Store at 4°C short term.Store at -20°C long term.Avoid freeze / thaw cycle.
Alias	DDDDK Tag, DYKDDDDK Tag
Image	<p>1ug Flag fusion protein+ Primary antibody dilution at 1:10000</p>  <p>IF analysis of 293 cells transfected with a Flag-tag protein,1:2000 dilution (blue DAPI,red anti-Flag)</p> 



PRODUCT USE LIMITATION

These products are intended for research use only.