

Specsure[®] Calf serum (CS)

Cat# ACE005

Store at -20°C.

INFORMATION

Product Name	Specsure [®] Calf serum (CS)
Cat. No.	ACE005
Size	1L
Product type	From Naturally Clotted or Defibrinated Blood
Bioburden	<2,000 cfu/mL
Endotoxin	Less than 5.0 EU/mL
pH	6.0 to 8.5
Osmoality	240 to 340 mOsm/Kg
Haemoglobin	Less than 25 mg/dL
Total Protein	Greater than 45 g/L
Retention Samples	In 10 mL LDPE vials

PRODUCT USE LIMITATION

These products are intended for research use only.