

dNTP Mix (PCR grade)

Storage at -20 °C

INTRODUCTION

The dNTP Mix is an equimolar mixture of high purity dATP, dTTP, dCTP, and dGTP, which can be directly used for PCR, reverse transcription (RT) and one-step RT-PCR.

CONTENTS

2.5 mM dNTP		10 mM dNTP	
EP1009	1 ml	EP1011	1 ml
EP1010	1ml x 5	EP1012	1ml x 5

QUALITY CONTROL

Purity > 99% (HPLC).

No exonuclease, endonuclease, and RNase activities.

PRODUCT USE LIMITATION

These products are intended for research use only.