

Tween 60

Cat# C2021

Storage at Room temperature

INFORMATION

Chemical Formula: $C_{64}H_{126}O_{26}$


Molecular Weight: 1311.68

CAS. No: 9005-67-8

Density: 1.044 g/mL at 25 °C

Alias: Polysorbate61, T-61, PEG-3 sorbitan stearate, PEG-6 sorbitan stearate, PEG-40 sorbitan stearate, PEG-60 Sorbitan stearate, polyoxyethylenesorbitan monostearate

CHEMICAL STRUCTURE



PRODUCT USE LIMITATION

These products are intended for research use only.