


# DTNBP1 Antibody

Cat# A05009PC

Upon receipt, store at -20°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	DTNBP1 Antibody
<b>Cat. No.</b>	A05009PC
<b>Size</b>	100ug
<b>Uniprot</b>	Q96EV8
<b>Product type</b>	Primary antibody
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Immunogen</b>	Recombinant protein of Human DTNBP1
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Tested applications</b>	ELISA, WB, IHC
<b>Application</b>	WB:1:500-1:2000,IHC:1:50-1:200
<b>Conjugation</b>	Non-conjugated
<b>Purification Method</b>	Affinity purification
<b>Storage buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Alias</b>	Dysbindin (Biogenesis of lysosome-related organelles complex 1 subunit 8) (BLOC-1 subunit 8) (Dysbindin-1) (Dystrobrevin-binding protein 1) (Hermansky-Pudlak syndrome 7 protein) (HPS7 protein), DTNBP1, BLOC1S8
<b>Image</b>	 <p>Western blot analysis of extracts of mouse brain cell lines, using DTNBP1 antibody.</p>

## PRODUCT USE LIMITATION

These products are intended for research use only.