



DUSP12 Antibody

Cat# A05010PC

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	DUSP12 Antibody
Cat. No.	A05010PC
Size	50 ul, 100 ul
Uniprot	Q9UNI6
Product type	Primary antibody
Species Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of Human DUSP12
Host	Rabbit
Clonality	Polyclonal
Tested applications	ELISA, WB, IHC
Application	ELISA:1:1000-1:5000, WB:1:200-1:1000, IHC:1:25-1:100
Conjugation	Non-conjugated
Purification Method	Antigen affinity purification
Storage buffer	PBS pH7.4, 0.05% NaN ₃ , 40% Glycerol.
Alias	dual specificity phosphatase 12
Image	 <p>Jurkat cells, Primary antibody: A05010PC (DUSP12 Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution</p>



The image is immunohistochemistry of paraffin-embedded Human brain tissue using A05010PC (DUSP12 Antibody) at dilution 1/25.

PRODUCT USE LIMITATION

These products are intended for research use only.