



Goat Anti-Rabbit IgG (H+L)- HRP

Cat# A1013S – 1 mg

Storage at -20 °C

INFORMATION

Product Name	Goat Anti-Rabbit IgG (H+L)- HRP
Cat. No.	A1013S
Product type	Secondary antibodies
Species Reactivity	Rabbit
Immunogen	Rabbit IgG, whole molecule
Host	Goat
Clonality	Polyclonal
Target	Rabbit IgG
Concentration	1 mg/ml
Tested applications	WB, ELISA, IHC
Application	WB : 1:4000-20000 ELISA : 1:4000-20000 IHC : 1:500-1:1000
Conjugation	HRP
Purification	The antibody was affinity-purified from goat antiserum by affinity-chromatography using epitope-specific immunogen.
Storage buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage instruction	Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

PRODUCT USE LIMITATION

These products are intended for research use only.