




# EPAS1 Antibody

Cat # A240012

Upon receipt, store at -20°C. Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	EPAS1 Antibody
<b>Cat. No.</b>	A240012
<b>Size</b>	100 ug
<b>Uniprot</b>	Q99814
<b>Product type</b>	Primary antibody
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Immunogen</b>	Synthesized peptide derived from Human EPAS-1 around the non-acetylation site of K385.
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Tested applications</b>	ELISA, WB
<b>Application</b>	WB:1:500-1:2000,ELISA:1:20000
<b>Conjugation</b>	Non-conjugated
<b>Purification Method</b>	Affinity purification
<b>Storage buffer</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Alias</b>	EPAS1; BHLHE73; HIF2A; MOP2; PASD2; Endothelial PAS domain-containing protein 1; EPAS-1; Basic-helix-loop-helix-PAS protein MOP2; Class E basic helix-loop-helix protein 73; bHLHe73;
<b>Image</b>	 <p>Western Blot analysis of Jurkat L929 cells using EPAS-1 Polyclonal Antibody</p>

## PRODUCT USE LIMITATION

These products are intended for research use only.

Technical support: [service@acebiolab.com](mailto:service@acebiolab.com)

Phone: 886-3-2870051

ACE BIOLABS CO., LTD. [www.acebiolab.com](http://www.acebiolab.com)

6F., No 68, Qingsheng 2nd St., Dayuan Dist., Taoyuan City33743, Taiwan.