



# Smad1/5/9 Polyclonal Antibody

Cat# A22916PI

Upon receipt, store at -20°C . Avoid repeated freeze.

## INFORMATION

Product Name	Smad1/5/9 Polyclonal Antibody
Cat. No.	A22916PI
Size	20 µl, 50 µl ,100 µl
Product type	Primary Antibody
Species Reactivity	Hu, Ms, Rt
Gene Id / Swiss Prot No	Human Gene Id : <a href="#">4086/4090/4093</a> Human Swiss Prot No : <a href="#">Q15797/Q99717/O15198</a> Mouse Gene Id : <a href="#">17125/17129/55994</a> Rat Gene Id : <a href="#">59328/85435</a> Rat Swiss Prot No : <a href="#">P97588/Q9R1V3/O54835</a>
Immunogen	The antiserum was produced against synthesized peptide derived from human Smad1/5/9. AA range: 10-59
Host	Rabbit
Clonality	Polyclonal
Tested applications	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. .
Conjugation	Non-Conjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration	1 mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Image	 <p>Western Blot analysis of HeLa cells using Smad1/5/9 Polyclonal Antibody diluted at 1 : 500</p>



### PRODUCT USE LIMITATION

These products are intended for research use only.