


Caspase-3 Polyclonal Antibody

Cat# A20591PI- 20/50/100 ug

Storage at -20 °C

INFORMATION

Product Name	Caspase-3 Polyclonal Antibody
Cat. No.	A20591PI
Product type	Primary antibodies
Species reactivity	Human, Mouse, Rat
Clonality	polyclonal
Host	Rabbit
Size	20 ug/50 ug/100 ug
Uniprot	Human: P42574/ Mouse: P70677/ Rat: P55213
Gene name	CASP3
Protein name	Caspase-3
Immunogen	The antiserum was produced against synthesized peptide derived from human Caspase 3. AA range:116-165
Specificity	Caspase-3 Polyclonal Antibody detects endogenous levels of Caspase-3 protein.
Concentration	1 mg/ml
Tested applications	WB, IHC-p, IF, ELISA
Application	WB : 1:500-1:2000/IHC-p : 1:100-1:300/IF : 1:200-1:1000/ELISA : 1:5000. Not yet tested in other applications.
Conjugation	None
Purification	Affinity purified
Molecular Weight	37 kDa
Stock buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Storage instruction	Store at -20°C or -80°C for 1 year. Avoid repeated freeze/thaw cycles.
Alias	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1



PRODUCT USE LIMITATION

These products are intended for research use only.