

# 5X RIPA Lysis Buffer

Cat# C1002 – 250 ml | C1003 – 500 ml

Store at 4 °C

## **INTRODUCTION**

RIPA Lysis Buffer is a 5X concentrated formulation reagent for cell lysis.

## **PRODUCT COMPONENTS**

Diluted to a 1X working concentration, contains	
Tris-HCl (pH 7.6)	25 mM
NaCl	150 mM
NP-40	1%
Sodium deoxycholate	1%
SDS	0.1%

## **RELATED PRODUCTS**

A1029 100X Protease Inhibitor Cocktail (EDTA-free)

A1030 100X Phosphatase Inhibitor Cocktail (A+B)

## **PRODUCT USE LIMITATION**

These products are intended for research use only.