

Potassium phosphate monobasic (KH₂PO₄)

Cat# C3112 – 500 g**Storage at Room Temperature****INFORMATION****CAS:** 7778-77-0**M.W.:** 136.09**Formula:** KH₂PO₄**Purity:** >99%**Appearance:** White crystal powder**Solubility:** Soluble in water, clear, colorless.**pH (50 g/L, 25°C):** 4.2~4.5

Trace: SO ₄ ²⁻	<0.01%
Cl ⁻	<0.005%
Pb	<0.003%
Fe	<0.005%
As	<3 ppm
N	<0.005%
Insoluble matter	<0.05%

PRODUCT USE LIMITATION

These products are intended for research use only.