

SDMA (Symmetric dimethylarginine)

Cat# C3188 - 5 mg | C3189 - 25 mg | C3190 - 100 mg

Storage at -20°C.

INFORMATION

Description: SDMA (Symmetric dimethylarginine) is an endogenous inhibitor of nitric oxide (NO) synthase activity

CAS: 30344-004

M.W.: 202.25


Formula: C₈H₁₈N₄O₂

Appearance: White to off-white powder

Purity: >98%

Solubility: Soluble in DMSO (50 mg/ml); Need ultrasonic and adjust pH to 3 with 1M HCl)

CHEMICAL STRUCTURE



PRODUCT USE LIMITATION

These products are intended for research use only.