


V5-Tag Monoclonal Antibody

Cat# A00030TC

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	V5-Tag Monoclonal Antibody
Cat. No.	A00030TC
Size	50ug, 100ug
Product type	Primary antibody
Species Reactivity	N/A
Immunogen	Synthetic Peptide
Host	Mouse
Clonality	Monoclonal
Tested applications	WB:1:500-1:5000
Application	ELISA, WB
Conjugation	Non-conjugated
Storage buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Image	 <p>Western blot analysis of V5 recombinant protein, diluted at 1:5000.</p>

PRODUCT USE LIMITATION

These products are intended for research use only.