

FOLR1 Antibody, Biotin conjugated

Cat# A04814PC

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	FOLR1 Antibody, Biotin conjugated
Cat. No.	A04814PC
Size	50ug, 100ug
Uniprot	P15328
Product type	Primary antibody
Species Reactivity	Human
Immunogen	Recombinant Human Folate receptor alpha protein (25-234AA)
Host	Rabbit
Clonality	Polyclonal
Application	ELISA
Conjugation	Biotin
Purification Method	>95%, Protein G purified
Storage buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Alias	Folate receptor alpha (FR-alpha), Adult folate-binding protein (FBP), Folate receptor 1, (Folate receptor, adult), KB cells FBP, Ovarian tumor-associated antigen MOv18, FOLR1, FOLR

PRODUCT USE LIMITATION

These products are intended for research use only.